

























Whisker Test Result

Package	Preconditioned	TEST	CONDITION	Maximum Whisker Measurement					
				1 (1 mth)	2 (2 mth)	3 (3 mth)	4 (4 mth)	5 (5 mth)	6 (6 mth)
16L SO	Preconditioned without IR reflow	Storage Test	20-30°C, 30-80%RH, 1-6months	0	0	0	0	0	0
									
		Aging Test	60°C, 85%RH, 1-6 months	0	0	0	0	0	0
									
	Preconditioned with IR reflow	Storage Test	20-30°C, 30-80%RH, 1-6months	0	0	0	0	0	0
									
		Aging Test	60°C, 93%RH, 1-6 months	0	0	0	0	0	0
									
		Temperature Cycling	-55°C to 85°C, 500cycles, 1000cycles	0	0				
				